


Bronchodilator Efficacy Delivered Using an In-line Aerogen® Solo during High-flow Nasal Cannula Therapy Versus a Jet Nebulizer and Facemask

Original article: Reminiac F, Vecellio L, Bodet-Contentin L, et al. Nasal high-flow bronchodilator nebulization: a randomized cross-over study. Ann Intensive Care. 2018;8(1):128.

Background



HFNC therapy is increasingly used as a method of non-invasive respiratory support; however, there is a lack of controlled clinical data on the efficacy of aerosolized bronchodilator therapy delivered in this setting

Objective




The aim of this study was to compare efficacy outcomes with an aerosolized bronchodilator delivered using the Aerogen Solo in-line with HFNC versus a standard jet nebulizer with facemask; for comparative purposes, outcomes were also assessed with HFNC therapy alone


Materials and Methods

Design: Randomized, crossover study

Adult patients with reversible obstructive lung disease (COPD/ asthma)*



Administered in a randomized order on 3 separate days across 1 week



Outcome measures




- The primary endpoint was the relative improvement in FEV₁ with the Aerogen Solo versus a jet nebulizer and facemask
- Bronchodilator efficacy was assessed using spirometry and plethysmography performed before and after each HFNC session

*FEV₁/FVC ratio <70% and a positive bronchodilator reversibility test (FEV₁ increase of ≥12% and ≥200 mL following administration of inhaled albuterol). FEV₁, forced expiratory volume in 1 second; FVC, forced vital capacity; HFNC, high-flow nasal cannula.

Bronchodilator Efficacy Delivered Using an In-line Aerogen® Solo Versus a Jet Nebulizer and Facemask During High-flow Nasal Cannula Therapy

Original article: Reminiac F, Vecellio L, Bodet-Contentin L, et al. Nasal high-flow bronchodilator nebulization: a randomized cross-over study. Ann Intensive Care. 2018;8(1):128.

Relative increase in FEV₁ following administration of aerosolized bronchodilator therapy



Changes in FEV₁ and residual volume with the Aerogen Solo and jet nebulizer were well correlated



Both delivery methods were well tolerated and had similar patient comfort ratings



Per the authors, the results demonstrate that the **Aerogen Solo** delivers inhaled therapy “within the [HFNC] circuit with **the same efficacy and tolerance** [as standard facemask jet nebulization] **avoiding cumbersome equipment switches**”



Bronchodilator efficacy was comparable with an in-line Aerogen Solo during HFNC therapy versus a jet nebulizer with facemask, with **no interruption of HFNC** therapy required with the Aerogen Solo

Want to know more?
Scan or click the QR code

