

Aerogen® Solo

Dry Side of the Humidifier

Quick Set-up Guide



Aerogen® Solo
AG-AS3200



T-Piece
AG-AS3010



Aerogen Pro-X Controller

Two modes of operation:

- 30 minutes
- Continuous




Aerogen USB Controller

Two modes of operation:

- 30 minutes
- 6 hours

Note: Humidifier parts not supplied by Aerogen.
Consult manufacturer for usage guidelines.

The Aerogen Solo delivers **4 times more medication** to patient's lungs.¹



✓ More efficient treatment of the patient²



1. ARI, A., ATALAY, O. T., HARWOOD, R., SHEARD, M. M., ALJAMHAN, E. A. & FINK, J. B. 2010. Influence of nebulizer type, position, and bias flow on aerosol drug delivery in simulated pediatric and adult lung models during mechanical ventilation. *Respir Care*, 55, 845-51.
2. Berlinski A, Willis JR. Albuterol delivery by 4 different nebulizers placed in 4 different positions in a pediatric ventilator in vitro model. *Respiratory care*. 2013;58:1124-1133