

Aerogen® Solo

Adult Nasal High Flow

Quick Set-up Guide



Aerogen® Solo
AG-AS3200



T-Piece
AG-AS3010



22F with
Oxygen Nipple



Aerogen Pro-X Controller

Two modes of operation:

- 30 minutes
- Continuous



Aerogen USB Controller

Two modes of operation:

- 30 minutes
- 6 hours

Note: Humidifier parts not supplied by Aerogen.
Consult manufacturer for usage guidelines.

- ✓ Effective drug delivery¹
- ✓ Less drug wastage²
- ✓ Better patient tolerance²



1. BHASHYAM, A. R., WOLF, M. T., MARCINKOWSKI, A. L., SAVILLE, A., THOMAS, K., CARCILLO, J. A. & CORCORAN, T. E. 2008. Aerosol delivery through nasal cannulas: an in vitro study. *J Aerosol Med Pulm Drug Deliv*, 21, 181-8.

2. REMINIAC, F., VECCELIO, L., HEUZE-VOURC'H, N., PETITCOLLIN, A., RESPAUD, R., CABRERA, M., LE PENNEC, D., DIOT, P. & EHRMANN, S. 2015. Aerosol Therapy in Adults Receiving High Flow Nasal Cannula Oxygen Therapy. *J Aerosol Med Pulm Drug Deliv*.