

Aerogen® Ultra

Spontaneous Breathing

Quick Set-up Guide



Aerogen Pro-X Controller

Two modes of operation:

- 30 minutes
- Continuous




Aerogen USB Controller

Two modes of operation:


- 30 minutes
- 6 hours

Note: Aerogen Ultra is qualified for up to 3 hours of continuous use.
Note: Humidifier parts not supplied by Aerogen.
Consult manufacturer for usage guidelines.

The Aerogen Ultra delivers
significantly more medication
in half the time.¹



Small Volume Nebuliser



Aerogen Ultra²

- ✓ Reduced treatment time¹
- ✓ More efficient treatment of the patient²
- ✓ Totally silent



1. Hickin S, Mac Loughlin R, Sweeney L, Tatham A, Gidwani S. Comparison of mesh nebuliser versus jet nebuliser in simulated adults with chronic obstructive pulmonary disease. Poster at the College of Emergency Medicine Clinical Excellence Conference. 2014

2. Alcoforado, L, et al, Analysis of Deposition Radioaerosol nebulizers membrane in Healthy Subjects. Poster at International Society for Aerosols in Medicine 2015